

# The Clarity RSV Antigen Test

The Clarity Diagnostics RSV Antigen test is designed for the qualitative detection of respiratory syncytial virus antigen fusion protein directly from nasopharyngeal swab or nasal aspirate specimens. The kit includes all components to run a nasopharyngeal specimen, no additional components necessary. Proudly made in the USA!

Product Code: CLA-RSV20



- ✓ All Necessary Components Included in Kit\*
- ✓ Approved for 6 & Under & Adults 60+
- ✓ CLIA Waived
- ✓ 20 Tests Per Box
- ✓ Results in 15 Minutes
- ✓ 24-month Shelf Life
- ✓ External Controls Included
- ✓ Made in the USA

**Kit Includes:** • 20 Individually Packaged Test Devices • 20 Individually Packaged Nasopharyngeal Swabs • 20 Extraction Tubes • 20 Specimen Processing Tubes • 20 Sterile Normal Saline Vials • 40 Disposable Pipettes • 1 x 2.0 mL R1 (Extraction Reagent Solution) • 1 x 1.0 mL PC (Positive Control Reagent) & • Instructions for Use.

\*For a nasopharyngeal specimen.



★ MADE IN ★  
**USA**

For more information please email  
[sales@ClarityDiagnostics.com](mailto:sales@ClarityDiagnostics.com) or call  
877-722-6339

[www.ClarityDiagnostics.com](http://www.ClarityDiagnostics.com)

**We Support Better  
Outcomes™**