

Clarity Drug Testing Dip Cards




- ✓ Easy-to-Use Dip Cards
- ✓ Multiple Drug Combinations Available
- ✓ 510K Cleared and CLIA Waived
- ✓ Results can be Photocopied



Drug of Abuse	Abbreviation	Cut-Off
Amphetamine	AMP	500 & 1000 ng/ml
Barbiturates	BAR	300 ng/ml
Buprenorphine/ Norbuprenorphine	BUP/NBUP	10 ng/ml
Benzodiazepine	BZO/BZD	300 ng/ml
Cocaine	COC	150 & 300 ng/ml
Methamphetamine	mAMP/MET	500 & 1000 ng/ml
MDMA or Ecstasy	MDMA	500 ng/ml
Morphine	MOP	300 ng/ml
Methadone	MTD	300 ng/ml
Opiate	OPI	2000 ng/ml
Oxycodone	OXY	100 ng/ml
Phencyclidine	PCP	25 ng/ml
Propoxyphene	PPX	300 ng/ml
Tricyclic Anti Depressants	TCA	1000 ng/ml
Marijuana/Hashish	THC	50 ng/ml

Read the Drug Test Results
at 5 Minutes.



*All positive results are presumptive and should be confirmed by an alternative method (e.g. GC/MS).



For more information please email
sales@claritydiagnostics.com or call 877-722-6339
www.claritydiagnostics.com

We Support
Better Outcomes™