

## Clarity Urine Liquid Controls Semiquantitative (Dipstick) Assay



- ✓ 22 Month Shelf Life
- ✓ Easy to Use Dropper Bottles
- ✓ Barcode Programmed with Lot Number & Expiration Date for Accuracy & Compliance

- ✓ Intended for use with Clarity Platinum Urinalysis 10 Parameter Strips (CLA-10P) and Clarity Visual Read Microalbumin/Creatinine Strips (CLA-MAC) Only
- ✓ Detects 10 Parameters present in Urine Sample including Microalbumin and Creatinine

### Contents:

- 1 X 25 ml Level1/Negative
- 1 X 25 ml Level2/Positive

Product Code:

**CLA-UHCRL25**



For more information please email  
[sales@claritydiagnostics.com](mailto:sales@claritydiagnostics.com) or call 877-722-6339

[www.claritydiagnostics.com](http://www.claritydiagnostics.com)

We Support  
Better Outcomes™