

Clarity Drug Testing Cups



- ✓ Built-In Temperature Strip
- ✓ Results in 5 Minutes
- ✓ Screw Top Lid for a No-Leak Cup
- ✓ Multiple Drug Combinations Available
- ✓ 510K Cleared and CLIA Waived



DRUG OF ABUSE	ABBREVIATION	CUTOFF
Amphetamine	AMP	500 ng/mL or 1000 ng/mL
Secobarbital	BAR	300 ng/mL
Secobarbital	BUP	10 ng/mL
Oxazepam	BZO	300 ng/mL
Cocaine	COC	150 ng/mL
2-Ethylidene-1,5-Dimethyl-3,3-Diphenylpyrrolidine	EDDP	300ng/ml
Fentanyl	FYL	1 ng/mL
Ecstasy	MDMA	500 ng/mL
Propoxyphene	PPX	300 ng/mL
Opiate	MOR/OPI	300 ng/ml or 2000 ng/ml
Methadone	MTD	300 ng/mL
Phencyclidine	PCP	25 ng/mL
Oxycodone	OXY	100 ng/mL
Norfentanyl	NFYL	5 ng/mL
Methamphetamine	MET	500 ng/mL or 1000 ng/mL
Nortriptyline	TCA	1000 ng/mL
6-Monoacetylmorphine	6-MAM	10 ng/mL
Tramadol	TML	100 ng/mL
Marijuana	THC	50 ng/mL



For more information please email sales@claritydiagnostics.com or call 877-722-6339
www.claritydiagnostics.com

We Support
 Better Outcomes™