

CLARITY

Together...We Support Better Outcomes™





TABLE OF CONTENTS




COVID Products

Our COVID-19 products are designed to deliver fast, easy-to-use, and highly accurate results.



URINALYSIS Products

We offer a complete set of Urinalysis products that are designed to provide you with quicker, reliable results.



INFECTIOUS DISEASE Products

We provide easy-to-use diagnostic products that can quickly test for infectious diseases like Mononucleosis and Strep A.



PREGNANCY Products

Clarity offers your practice the opportunity to bring your pregnancy testing in-house with our array of products.



GLUCOSE MONITORING Products

For patients with diabetes, glucose monitoring is essential to ensure their well being.



DRUGS OF ABUSE Products

Allows your facility to quickly and easily test patients and identify accurately elevated levels of commonly abused drugs.



HEMOGLOBIN Meters

Measuring hemoglobin is essential for checking the health of your patients and maintaining their quality of life.

CLARITY

PLATINUM URINE ANALYZER

Discover the benefits of our Clarity Platinum analyzer. Proudly made in the USA, the analyzer comes pre-configured for plug-and-play operation with a color touch screen and user-friendly interface.




HIGH QUALITY COLOR
TOUCH SCREEN
USER FRIENDLY
INTERFACE



DIRECT CLOUD
CONNECTIVITY

CLA-PLTUA



- ✓ INCLUDES NEW STRIP LOT DETECTION TECHNOLOGY
- ✓ ADDITIONAL SENSOR TO DETECT HUMIDITY OVEREXPOSURE

- ✓ AUTOMATIC CALIBRATION
- ✓ BUILT-IN THERMAL PRINTER WITH AUTO FEEDER

- ✓ STORES 2000 TEST RESULTS
- ✓ TECHNOLOGY TO PREVENT USE OF EXPIRED STRIPS

- ✓ BARCODE SCANNER INCLUDED

CLARITY PLATINUM MARKET ANALYSIS AND COMPETITIVE ADVANTAGES

	 CLARITY PLATINUM	 SIEMENS CLINTEK STATUS	 MISSION D120 SMART	 BIOSYS URISCAN	 HENRY SCHEIN URISPEC PLUS	 MCKESSON CONSULT 120	 ROCHE URISYS 1100
Printer w/ Autofeeder	✓	X	X	X	X	X	X
Results Memory	2,000	700	2,000	2,000	200	2,000	100
LCD Display	Color LCD Touchscreen	Single Color Touchscreen	Single Color Touchscreen	Color LCD Touchscreen	Single Color Touchscreen	Single Color Touchscreen	Single Color Touchscreen
Operator Login	✓ to 500 operators	Up to 700 operators	Up to 800 operators	X	X	X	X
Admin Login	✓ Up to 100 admins	X	X	X	X	X	X
Compliant Strip Design	✓ IR Band	Auto Check Pad	X	X	X	X	✓
Humidity Check	✓	✓	X	X	X	X	X
Strip Expiration Compliance	✓ Via Bar code Scan	X	X	X	X	X	X
Control Expiration Compliance	✓	X	X	X	X	X	X
Wrong Control Level Warning	✓	X	X	X	X	X	X
Direct Cloud Connectivity	✓	X - Additional Cost	X	X	X	X	X
Data Output	RJ45 Ethernet	RS-232	RS-232	RS-232	RS-232 / RJ45 (Inactive)	RS-232	RS-232
Barcode Reader Included	✓	X	X	X	X	X	X
	✓	✓	X	X	X	X	X

CLARITY

PLATINUM 10 PARAMETER URINALYSIS STRIPS

Our CLA-10P urinalysis strips are for use only with the Clarity Platinum and are more than 99% accurate, and test for Leukocytes, Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ketone, Bilirubin, and Glucose.

Each strip includes a IR band to ensure compliance and our lid-lock technology minimizes the strips' exposure to air and humidity.



**MORE THAN 99%
ACCURATE WITH
CONSISTENT RESULTS**

CLA-10P



LEUKOCYTES	NEGATIVE 0	TRACE 15	SMALL (+) 70	MODERATE (++) 125	LARGE (+++) 500 ca CELLS/ μ L		
NITRITE	NEGATIVE	POSITIVE (ANY DEGREE OF PINK COLOR)					
URO-BILINOGEN	0.2 NORMAL 3.2	1 16	2 32	4 64	8 128	EHRlich UNITS/dl URINE μ mol/L	
PROTEIN	NEGATIVE	TRACE	30 (+) 0.3	100 (++) 1.0	300 (+++) 3.0	mg/dl g/L	
PH	5.0	6.0	6.57	7.0	7.58	8.0	8.5
BLOOD	NEGATIVE 0	NON-HEMOLYZED 10 TRACE	HEMOLYZED TRACE	SMALL (+) 25	MODERATE (++) 80	LARGE (+++) 200 ca CELLS/ μ L	
SPECIFIC GRAVITY	1.000	1.005	1.010	1.015	1.020	1.025	1.030
KETONE	NEGATIVE	TRACE (5) 0.5	SMALL (15) 1.5	MODERATE (40) 4.0	LARGE (80) 8.0	mg/dl mmol/L	
BILIRUBIN	NEGATIVE	SMALL +	MODERATE ++	LARGE +++			
GLUCOSE	NEGATIVE	100 5	250 (+) 15	500 (++) 30	1000 (+++) 60	mg/dl mmol/L	



- ✓ FOR USE WITH CLARITY PLATINUM URINE ANALYZER ONLY
- ✓ BUILT-IN TECHNOLOGY PREVENTS USE OF EXPIRED STRIPS

- ✓ TESTS FOR LEUKOCYTES, NITRITE, UROBILINOGEN, PROTEIN, PH, BLOOD, SPECIFIC GRAVITY, KETONE, BILIRUBIN, AND GLUCOSE

- ✓ STRIPS WITH IR BAND TO ENSURE COMPLIANCE

CLARITY

PLATINUM URINE LIQUID CONTROLS

Our Clarity Platinum quality control sets are designed for use with the Clarity Platinum urine analyzer. They feature easy to use 25 mL dropper bottles, and their labels include barcodes with lot and expiry data for accuracy and compliance.



CLA-UHCRL25



 BARCODE PROGRAMMED FOR ACCURACY AND COMPLIANCE

CLARITY



COVID-19 ANTIGEN RAPID TEST 20 PACK

In response to the worldwide pandemic, Clarity is proud to introduce a new line of reliable COVID-19 products.

Our COVID-19 products are proudly manufactured in the USA. With a continuing pandemic, the need for accurate results and reliable tests are of the utmost importance, which is why our COVID products are made to be quick, easy to use, and above all accurate.



CLA-COV-AG20



✓ 510K CLEARED

✓ POCT APPROVED

✓ RELATIVE SPECIFICITY: >99.69%

✓ CLIA WAIVED

✓ ANTERIOR NASAL SWAB TEST

✓ RESULTS IN 15 MINUTES



RESULTS IN
15 MINUTES

CLARITY

UROCHECK URINE STRIPS

Clarity Urocheck urine reagent strips are for visual read use only. The Urocheck line includes a wide variety of models for different clinical needs, including 2GP, 5OB, and even microalbumin and creatinine strips.

DTG-10SG | DTG-5GP
DTG-105G | DTG-50B



100 STRIPS
PER BOX



✓ CONSISTENT RESULTS

✓ CAP/API VALUES LISTED

✓ ROOM TEMPERATURE STORAGE

✓ TESTS IN THE SAME ORDER
AS LEADING BRANDS

CLARITY

CD-UCTL30



UROCHECK URINE LIQUID CONTROLS


The Clarity CD-UCTL30 Liquid Controls are for use with the Clarity Urocheck brand urine reagent strips. They ship in 3 sets of 10 mL easy-to-use dropper bottles with color coded caps.



EASY DISPENSING
DROPPER BOTTLES



22 MONTH
SHELF LIFE



THEY DETECT
ALL 10 PARAMETERS
OF OUR DTG-10SG



IDEAL FOR POC
APPLICATIONS



✓ **STABLE IN A LIQUID STATE FOR 2 YEARS, EVEN AFTER BEING OPEN!**

✓ **VALUES FOR ALL ANALYTES, INCLUDING CREATININE AND MICROALBUMIN, ARE PROVIDED**

✓ **INTENDED FOR USE WITH UROCHECK URINALYSIS STRIPS**

CLARITY

TOX ALERT FIA READER

The Clarity Tox Alert FIA reader is an incredibly sensitive fentanyl testing device that detects the presence of fentanyl in urine to a cutoff of 1ng/mL.

It is CLIA waived and OTC approved, and includes a wireless printer via Bluetooth connection for secure and seamless printing of test results.



RESULTS IN
5 MINUTES



EASY TO USE
MINIMAL TRAINING
REQUIRED



HANDHELD:
LIGHT WEIGHT (5.76 OZ)



✓ ULTRA-SENSITIVE: 1 NG/ML

✓ OVERALL ACCURACY: >99%

✓ PRINTABLE RESULTS VIA BLUETOOTH PRINTER


✓ SMALL SAMPLE (3 DROPS) OF URINE

✓ REIMBURSABLE WITH ANALYZER
READ CPT CODE-80306-QW



CLARITY

DTG-STP25



STREP A RAPID TEST STRIPS

Clarity's DTG-STP25 Strep A rapid test kit is made in the USA and has a sensitivity of 95% and a specificity of 98.4%.

DTG-STP25 is CLIA waived and FDA cleared, with controls included as part of the kit.



✓ SENSITIVITY OF 95%


✓ SPECIFICITY OF 98.4%

✓ ROOM TEMPERATURE STORAGE

✓ CLIA WAIVED

✓ STRIPS IN VIAL FORMAT

✓ CONTROLS INCLUDED IN KIT



RSV ANTIGEN TEST

20 PACK

Clarity's CLA-RSVN525 RSV test is made in the USA and provides rapid results within 15 minutes. The kit includes everything needed to run the test, including controls.

The sample requires an anterior nasal swab for easier collection. For use with children between 6 months and 6 years old and adults that are 60 years old and beyond.



✓ ALL NECESSARY COMPONENTS INCLUDED IN KIT*


✓ EXTERNAL CONTROLS INCLUDED
✓ ANTERIOR NASAL SWAB TEST

✓ APPROVED FOR 6 & UNDER & ADULTS 60+

CLARITY

MONONUCLEOSIS TEST KIT

Clarity's Mononucleosis rapid test kit has 99% overall accuracy and is made in the USA. It can be used with both whole blood (CLIA waived) and plasma/serum (CLIA moderate).



DTG-MONO



✓ SENSITIVITY OF 99%

✓ SPECIFICITY OF 98.6%

✓ CLIA WAIVED FOR WHOLE BLOOD ONLY,
CAN BE USED WITH PLASMA/SERUM

✓ OVERALL ACCURACY: >99%

✓ SERUM IS CLIA MODERATE

CLARITY

DTG-PLUS25



HCG SINGLE STEP URINE CASSETTE PREGNANCY TEST

25 **PACK**

Clarity's DTG-PLUS25 pregnancy test cassette has both 99% overall accuracy and specificity of 100%. Results within as little as 5 minutes. Each box includes 25 cassettes and a product insert.

The DTG-UHCTLS hCG controls are sold separately.



18 MONTH DATING

CLIA
WAVED



SPECIFICITY OF 100%
NO CROSS REACTION
WITH FSH/TSH/LH



SENSITIVITY OF 25 mIU/ml



OVERALL ACCURACY: >99%



CLIA WAIVED

CLARITY

DTG-COMBO



PREGNANCY COMBO (URINE/SERUM) TEST CASSETTES


Clarity's DTG-COMBO pregnancy test cassette combines urine and serum detection of hCG into one product. The test has an overall accuracy of 99% with specificity of 100%, at a cutoff of 20 mIU/mL for urine and 10 mIU/mL for serum. Urine is CLIA waived, serum is CLIA moderate. Uses the DTG UHCTLS hCG controls, which are sold separately. Each box includes 50 tests.



 24 MONTH
DATING

 50
50 TESTS PER BOX
+ PACKAGE INSERT

 BUILT-IN QUALITY CONTROL
 OVERALL ACCURACY: >99%

 SENSITIVITY:
10 U/ml - Serum
20 U/ml - Serum

 ROOM TEMPERATURE STORAGE

 SPECIFICITY OF 100%

CLARITY

DTG-HCG25



PREGNANCY TEST STRIPS

Clarity's hCG pregnancy dipstick test kits are one-step pregnancy test strips with 25 mIU/mL sensitivity with 99% overall accuracy and 100% specificity. Each box includes 25 test strips and a package insert.

The DTG-UHCTLS hCG controls for running quality control are sold separately



✓ ONE-STEP TEST

✓ OVERALL ACCURACY: >99%

✓ ROOM TEMPERATURE STORAGE

✓ BUILT-IN QUALIT CONTROL

✓ SENSITIVITY: 25 mIU/ml

✓ SPECIFICITY OF 100%

CLARITY

URINE/HCG LIQUID CONTROLS SEMIQUANTITATIVE (DIPSTICK) ASSAY

Clarity's DTG-UHCTLS hCG and urine controls can be used for quality control with many of the most popular brands in urinalysis as well as hCG pregnancy tests.

The easy-to-use dropper bottles allow precise and controlled administration of the solution to save on messes and wasted product.



DTG-UHCTLS



INTENDED FOR USE WITH:

- SIEMENS®,
- MULTISTIX® STRIPS,
- CLINITEK® SERIES OF ANALYZERS,
- CHEMSTRIP® DIPSTICKS,
- CHEMSTRIP® MICRAL,
- ROCHE URISYS 1100 AND MINI UA URINE ANALYZER SERIES*.


AS WELL AS ALL CLARITY DIAGNOSTICS HCG PREGNANCY TESTS, AND WITH OTHER COMMERCIALY AVAILABLE EARLY PREGNANCY TESTS.

✓ INCLUDES HCG AND CAN BE USED FOR QUALITATIVE RESULTS IN RAPID EARLY DETECTION HCG (PREGNANCY) TESTS

✓ INCLUDES VALUES FOR ALL ANALYTES, CREATININE, MICROALBUMIN AND HCG.

CLARITY

BLOOD GLUCOSE MONITORING SYSTEM



CD-BG4 CONTENTS:

- ✓ 1 EACH OF NORMAL AND HIGH CONTROL
- ✓ PACKAGE INSERT

DTG-GL6 CONTENTS:

- ✓ 1 LANCET DEVICE

DTG-GL7 CONTENTS:

- ✓ SINGLE USE TWIST TOP STERILE 30G LANCETS, 100/BX

DTG-SL23100 CONTENTS:

- ✓ STERILE SINGLE USE 26G SAFETY LANCETS, 100/BX

CD-BG5 CONTENTS:

- ✓ 50 TEST STRIPS
- ✓ PACKAGE INSERT

CLARITY

CD-CDOA-8145FNF
CD-CDOA-8144FNF

CD-CDOA-6125
CD-CDOA-1144

CD-CDOA-4104
CD-CDOA-7125
CD-CDOA-6124



DRUG TESTING CUPS

Clarity's Drugs of Abuse cups include a temperature strip, are 510K cleared and CLIA waived, and are available in multiple combinations, from 5 to as many as 19 panels. Results can be interpreted in as little as 5 minutes, and the included screw top lid prevents leakage. We can create custom panels and offer Forensic Use Only options as well.



Drug of abuse	Abbreviation	Cut off
Amphetamine	AMP	500 ng/mL or 1000 ng/mL
Secobarbital	BAR	300 ng/mL
Buprenorphine	BUP	10 ng/mL
Oxazepam	BZO	300 ng/mL
Cocaine	COC	150 ng/ml or 300 ng/ml
2-Ethylidene-1,5-Dimethyl-3,3-Diphenylpyrrolidine	EDDP	300 ng/mL
Fentanyl	FYL	1 ng/mL
Ecstasy	MDMA	500 ng/mL
Propoxyphene	PPX	300 ng/mL
Opiate	MOR/OPI	300 ng/mL or 2000 ng/mL
Methadone	MTD	300 ng/mL
Phencyclidine	PCP	25 ng/mL
Oxycodone	OXY	100 ng/mL
Norfentanyl	NFYL	5 ng/mL
Methamphetamine	MET	500 ng/mL or 1000 ng/mL
Nortriptyline	TCA	1000 ng/mL
6-Monoacetylmorphine	6-MAM	10 ng/mL
Tramadol	TML	100 ng/mL
Marijuana	THC	50 ng/mL

BUILT-IN TEMPERATURE STRIP

MULTIPLE DRUG COMBINATIONS AVAILABLE

SCREW TOP LID FOR A NO-LEAK CUP

CLARITY

FUO DRUG TESTING CUPS

Clarity's Drugs of abuse cups include a temperature strip and are available in multiple combinations, from 5 to as many as 29 panels. Results can be interpreted in as little as 5 minutes and the included screw top lid prevents leakage.



SAMPLE DRUG COMBINATION:

Drug	Calibrator	Cut-Off (ng/mL)
Ethyl Glucuronide (EtG)	Ethyl Glucuronide	300
Fentanyl (FYL)	Fentanyl	20
Gabapentin (GAB)	Gabapentin	1,000
Synthetic Cannabinoid (K2)	JWH-018 5-pentanoic acid JWH-073 4-butanoic acid	50
Ketamine (KET)	Ketamine	1,000
Kratom (KRA)	Mitragynine	500
Psilocybin (PY)	Psilocin	500
Marijuana (THC)	11-nor- Δ^9 -THC-9 COOH	25 / 50 / 200
Tianeptine (TIA/ZAZA)	Tianeptine	500
Xylazine (XYL)	Xylazine	1,000
Adulteration Strip	Creatinine, Specific Gravity, pH	—

 BUILT-IN TEMPERATURE STRIP

 MULTIPLE DRUG COMBINATIONS AVAILABLE

 SCREW TOP LID FOR A NO-LEAK CUP



CLARITY

CD-DOA-1115
CD-DOA-6125

CD-DOA-154
CD-DOA-254
CD-DOA-1104



DRUG TESTING DIP CARDS

Clarity's multidrug dip cards are easy-to-use, 510K cleared and CLIA waived, and include all of the most popular combinations to fit any kind of use-case. We can create custom panels and offer Forensic Use Only options as well.



TEST TIME:
5 MINS



✓ CLARITY DRUG TESTING DIP CARDS

✓ EASY-TO-USE DIP CARDS

✓ RESULTS CAN BE PHOTOCOPIED

✓ CUSTOM MULTI-DRUG CUPS AND PANELS AVAILABLE WITH SPECIALTY DRUGS

✓ MULTIPLE DRUG COMBINATIONS AVAILABLE

✓ OVER 20 ADDITIONAL PARAMETERS AVAILABLE FOR FORENSIC USE ONLY

CLARITY

CD-UCTS90



STERILE URINE CONTAINER

The Clarity Sterile Urine Container is a 90 mL specimen cup with clear graduations in ounces and milliliters, a secure dual click cap for 100% leak-proof sealing, and a patient ID label included. Pneumatic tube safe up to 95 kPa and conveniently packed 400 per case.



400 CUPS
PER CASE



100%
LEAK-PROOF

✓ PATIENT ID LABEL INCLUDED

✓ TWIST-TOP DUAL CLICK CAP





✓ GRADUATIONS IN BOTH OUNCES AND MILLILITERS

✓ 90 ML CONTAINER WITH GRADUATIONS

✓ PNEUMATIC TUBE SAFE WITH 95KPA
SPECIFICATION

CLARITY

SALIVA DRUG TEST



E&I EXEMPT PANELS

Drug Tested	Cut-off
Alcohol (ALC)	0.02 g/dL
Amphetamine (AMP)	50 ng/mL
Barbiturates (BAR)	50 ng/mL
Benzodiazepines (BZO)	20/50 ng/mL
Buprenorphine (BUP)	5/10 ng/mL
Oxycodone (OXY)	20 ng/mL
Cocaine (COC)	20 ng/mL
Fentanyl (FTY)	20 ng/mL
Ketamine (KET)	50/100 ng/mL
Methadone (MTD)	30 ng/mL
Methamphetamine (MET/mAMP)	50 ng/mL
Methylenedioxy methamphetamine (MDMA)	50/100 ng/mL
Opiate (OPI)	40 ng/mL
Phencyclidine (PCP)	10 ng/mL
Synthetic cannabinoid (K2)	20 ng/mL
Marijuana (THC)	40 ng/mL

FORENSIC USE PANELS

Drug Tested	Cut-off
Amphetamine (AMP)	50 ng/mL
Barbiturates (BAR)	50 ng/mL
Cocaine (COC)	20 ng/mL
Marijuana (THC)	40 ng/mL
Methadone (MTD)	30 ng/mL
Methamphetamine (MET)	50 ng/mL
Methylenedioxy methamphetamine (MDMA)	50/100 ng/mL
Opiate (OPI)	40 ng/mL
Oxycodone (OXY)	20 ng/mL
Phencyclidine (PCP)	10 ng/mL



- ✓ ALL-IN-ONE MIDSTREAM DESIGN, SCREENS UP TO 14 DRUGS INCLUDING ALCOHOL
- ✓ EASY TO USE, NO ADDITIONAL EQUIPMENT NEEDED
- ✓ SIMPLE, OBSERVED COLLECTION PROCESS

- ✓ FORENSIC-USE-ONLY
- ✓ EMPLOYMENT AND INSURANCE EXEMPT PANELS AVAILABLE
- ✓ DETECTS RECENT DRUG USE

CLARITY

CD-DOACTLS10



DRUGS OF ABUSE CONTROLS

Clarity's Drugs of Abuse controls are made in the USA and are ready to use. They can be used for all Clarity dip and cup products, as well as other manufacturers' DOA products.

Each set includes a 10 mL negative control vial, a 10 mL positive control vial, and a package insert.



18 MONTH SHELF LIFE



- ✓ CONTROL MATRIX BASED ON HUMAN URINE TO SIMULATE PATIENT SAMPLES
- ✓ READY-TO-USE CONTROLS
- ✓ IDEAL FOR QUALITY CONTROL TESTING AND TRAINING

- ✓ POSITIVE CONTROL TESTS FOR MULTIPLE OPIATE CUTOFF LEVELS
- ✓ 510K CLEARED AND DEA EXEMPT
- ✓ SPECIALLY PACKAGED TO RETAIN POTENCY OF DRUGS

- ✓ INTENDED FOR USE WITH ANY BRAND OF DRUG CUP OR DIP CARD, AND ANY CLARITY DOA PRODUCT

CLARITY

HB CHECK HEMOGLOBIN TESTING SYSTEM

Our Hemoglobin Hb Check meter is CLIA waived, includes a large visual display, provides results in 15 seconds, and features a portable, handheld size. It tests both Hemoglobin and estimated Hematocrit. The unit requires minimal maintenance, and uses finger-stick whole blood samples. This device is not for home use/OTC use.



CLA-HB2



✓ FINGER-STICK WHOLE BLOOD SAMPLES

✓ APPROVED FOR MULTI-PATIENT USE

✓ MINIMAL MAINTENANCE

✓ SMALL, 1 μ L SAMPLE SIZE

✓ FOR TESTING HEMOGLOBIN &
ESTIMATED HEMATOCRIT IN BLOOD

✓ THIS DEVICE IS NOT FOR HOME USE


✓ 10 TEST CARTRIDGES PER VIAL

✓ 90 DAYS OPEN VIAL
STABILITY OF TEST CARTRIDGES



COMPLETE 12 LEAD PROFESSIONAL GRADE EKG SYSTEM FOR OPTIMAL PATIENT CARE

Our advanced EKG system is designed to ensure the highest standard of patient care, offering precision, ease of use, stickerless lead technology, and comprehensive data analysis that seamlessly integrates with your EMR.



✓ 12-LEAD STRAP SYSTEM

✓ NO ADHESIVE ELECTRODES

✓ CONSISTENT, ACCURATE PLACEMENT

✓ ADJUSTABLE FIT FOR ALL PATIENTS

✓ ANTIBACTERIAL & HYPOALLERGENIC

✓ COMPACT, EASY TO TRANSPORT

✓ RAPID DEPLOYMENT IN ANY CARE SETTING

✓ FAST AND SIMPLE OPERATION

✓ 80% FASTER THAN A TRADITIONAL EKG

✓ ONE-TOUCH RECORDING

✓ INTUITIVE APP

✓ SECURE CLOUD & EMR INTEGRATION

✓ QUICK HISTORICAL COMPARISON

✓ REAL-TIME DATA SHARING


CLARITY

FOLEY-SAFE 2.0

A catheter stabilizer for prevention of Catheter-Associated Urinary Tract Infections (CAUTI), reduction of pain, and prevention of unintended and accidental Foley catheter extractions.

Reduces pain and improves self-confidence and the patient's quality of life.

Foley-Safe 2.0 innovative technology guarantees that upon unintended or accidental pulling of the catheter, the catheter is locked within the device, or in the case of very strong pulling the catheter is severed, preventing transmission of these forces to the urethra, bladder neck, and the bladder itself.



Foley-Safe 2.0 provides three protection layers for patients with indwelling urethral catheters:

- ✓ **FIXATION:** The adhesive patch secures the catheter in place, preventing the urine bag from pulling the catheter.
- ✓ **ANTI-SLIDE SMART MECHANISM:** Detects and absorbs minor catheter movements and prevents the catheter from being pulled, thus potentially decreasing the chance of the catheter associated urinary tract infection (CAUTI).
- ✓ **DISCONNECTION:** When the pulling force is significant and dramatic, the advanced cutting mechanism is engaged smoothly in severing the catheter.