

Clarity Strep A

Rapid Test Strips

Product Code: DTG-STP25

Clarity Strep A is a rapid immunochromatographic assay for the qualitative detection of Group A Streptococcal antigen from throat swab specimens, aiding in the diagnosis of Group A strep infections.



Contents of the Box

- Strep A Test Strips in a Vial
- Disposable Extraction Tubes
- Reagents A & B
- Positive Control
- Negative Control
- Workstation
- Procedure Card
- Package Insert



Sensitivity: 95%



Specificity: 98.4%



24 Month Dating Open Vial Stability



Room Temperature Storage



Controls Included in Kit



Strips in Vial Format



Results in as fast as 5 minutes



Made in the USA



For more information, please email
sales@claritydiagnostics.com
or call (877) 722-6339



www.ClarityDiagnostics.com

