

Clarity Mono Cassettes



- | | |
|---|---|
| <ul style="list-style-type: none"> ✓ Sensitivity >99% ✓ Specificity 98.6% ✓ 99% Overall Accuracy ✓ Results as Early as 5 Minutes ✓ No Age Restriction | <ul style="list-style-type: none"> ✓ 18 Month Dating ✓ Room Temperature Storage ✓ Test Can Run with Capillary Finger-Stick Whole Blood, Plasma, or Serum |
|---|---|



Product Code:
DTG-MONO (15/bx)

Contents:

- Mono Test Cassettes
- Developer Solution
- Sample Transfer Tubes
- Positive Control
- Negative Control
- Instructional Insert
- Procedure Card



For more information please email
sales@claritydiagnostics.com or call 877-722-6339
www.claritydiagnostics.com

**We Support
Better Outcomes™**